

EQP-30

MANUAL HYDRAULIC PELLET PRESS



Getting reliable and repeatable results in RX or LIBS depends greatly on how the pellet is manufactured.

The EQP-30 Manual Hydraulic Pellet Press has been especially designed to prepare high quality pellets, homogeneous and stable.

This press can exert a pressure of up to 30 tons. It is a durable unit, safe and user friendly.



- **Robust and durable**
- **Prepare pellets of different diameters**
- **Maximum recommended pressure 30 tons**
- **Automatic release if excess pressure**

Technical specifications

Method:	pressure
Applications:	cement industry, metallurgy, thermal plants, environmental laboratories, materials recovery facilities, recycling industry, geology and mineralogy, ceramics
Maximum pressure:	30 tons
Dimensions:	410mm (high) x 315mm (long) x 200mm (wide)
Approx weight:	56Kg

© 2021 Equilab S.A.
All Rights Reserved

EQUILAB is a registered
trademark of Equilab S.A.



Equilab S.A.

Avda. Camino de lo Cortao, 21 - Nave 6
28703 San Sebastián de los Reyes - Madrid
Tel: 91 661 00 22 / Fax: 91 661 81 46
www.equilab.es
Customer care: equilab@equilab.es